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Adherence to recommended clearing, sterilization and lubrication practices will help to ensure the best product performance.

Important Suggestion on instruments' maintenance and sterilization

- 1) Before using the instrument they have to be washed with fresh water using only soft bristle toothbrush.
- 2) All the instruments presenting screws, joints or gains, such as extracting forceps, scissors needle hold etc., must always be open before washing and disinfecting.
- 3) Taking into consideration the fact that soap and disinfectant residuals could remain stuck on the instruments, rinse them off with demineralised or distilled water.
- 4) Immediately after this kind of operations, and before starting the sterilization process, dry the instruments very carefully.
- 5) All the instruments presenting screws, joint or gains, have to be carefully oiled using surgical instruments dedicated products or 100% pure "Paraffin oil" (never use industrial lubricants).

Sterilization process:

- Follow the temperatures suggested by the Autoclave producer.
- Put one instrument in each case before putting them in the Autoclave.
- Always open the instruments presenting joint or gains before starting sterilization.

"COLD PROCESSING" solutions:

- The effectiveness of the kind of solution must be carefully evaluated because they do not guarantee 100% sterilization against many bacteria and viruses.
- Taking into consideration the fact that sterilization needs many hours to be effective, the immersion in the solution could cause damage to the instruments and jeopardize their future use. We suggest, therefore, to reduce immersion time.

FV2000 Surgical Instruments is not responsible for damages caused to the instruments and due to an erroneous sterilization process.